



www.emrsafetyandhealth.com  
972-235-8330

# MEGACODE MEMORY

PULSE	NO PULSE	PCAC												
Vomit (Vitals, O2, Monitor, IV, Team)	Comit (CPR, O2, Monitor, IV, Team)	Post Cardiac Arrest Care												
<b>Bradycardia</b>	<b>Asystole/PEA</b>	<b>Vomit</b>												
Treat 2X Atropine Epi Pace	CPR Epi 3-5 min	(Vitals, O2, Monitor, IV, Team) • Adv Airway • 12 lead												
<b>Tachycardia</b>	<b>VF/VT</b>	<b>Targeted Temp Management</b>												
<table border="0"> <tr> <td style="vertical-align: top;"> <b>Stable</b>  <ul style="list-style-type: none"> <li>• Vagal</li> <li>• Adeno</li> <li>• Adeno</li> <li>• Sync</li> </ul> </td> <td style="vertical-align: top;"> <b>Unstable</b>  <ul style="list-style-type: none"> <li>• Sync</li> </ul> </td> </tr> </table>	<b>Stable</b> <ul style="list-style-type: none"> <li>• Vagal</li> <li>• Adeno</li> <li>• Adeno</li> <li>• Sync</li> </ul>	<b>Unstable</b> <ul style="list-style-type: none"> <li>• Sync</li> </ul>	32° - 37.5°C (36 HRS)	<b>H's T'S</b>										
<b>Stable</b> <ul style="list-style-type: none"> <li>• Vagal</li> <li>• Adeno</li> <li>• Adeno</li> <li>• Sync</li> </ul>	<b>Unstable</b> <ul style="list-style-type: none"> <li>• Sync</li> </ul>													
<b>High Quality CPR</b>	<table border="0"> <tr> <td style="border: 1px solid red; padding: 5px;"> <b>Shock</b> CPR 2-min Check Pulse VF/VT?         </td> <td style="border: 1px solid red; padding: 5px;">           VF/VT?  <b>Shock</b> CPR 2-min <b>Ami</b> Check Pulse         </td> </tr> <tr> <td style="border: 1px solid red; padding: 5px;"> <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse VF/VT?         </td> <td style="border: 1px solid red; padding: 5px;">           VF/VT?  <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse         </td> </tr> <tr> <td style="border: 1px solid red; padding: 5px;"> <b>Shock</b> CPR 2-min <b>Ami</b> Check Pulse VF/VT?         </td> <td style="border: 1px solid red; padding: 5px;">           VF/VT?  <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse         </td> </tr> <tr> <td style="border: 1px solid red; padding: 5px;"> <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse         </td> <td style="border: 1px solid red; padding: 5px;">           VF/VT?  <b>Shock</b> CPR 2-min <b>Nothing</b> Check Pulse         </td> </tr> </table>	<b>Shock</b> CPR 2-min Check Pulse VF/VT?	VF/VT? <b>Shock</b> CPR 2-min <b>Ami</b> Check Pulse	<b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse VF/VT?	VF/VT? <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse	<b>Shock</b> CPR 2-min <b>Ami</b> Check Pulse VF/VT?	VF/VT? <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse	<b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse	VF/VT? <b>Shock</b> CPR 2-min <b>Nothing</b> Check Pulse	<table border="0"> <tr> <td style="vertical-align: top; font-size: 2em;"><b>H's</b></td> <td style="border: 1px solid red; padding: 10px; text-align: center;">           Hypovolemia Hypoxia Hydrogen ion (acidosis) Hypo/Hyperkalemia Hypothermia         </td> </tr> <tr> <td style="vertical-align: top; font-size: 2em;"><b>T's</b></td> <td style="border: 1px solid red; padding: 10px; text-align: center;">           Tension Pneumothorax Tamponade (cardiac) Toxins Thrombosis (pulmonary) Thrombosis (coronary)         </td> </tr> </table>	<b>H's</b>	Hypovolemia Hypoxia Hydrogen ion (acidosis) Hypo/Hyperkalemia Hypothermia	<b>T's</b>	Tension Pneumothorax Tamponade (cardiac) Toxins Thrombosis (pulmonary) Thrombosis (coronary)
<b>Shock</b> CPR 2-min Check Pulse VF/VT?	VF/VT? <b>Shock</b> CPR 2-min <b>Ami</b> Check Pulse													
<b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse VF/VT?	VF/VT? <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse													
<b>Shock</b> CPR 2-min <b>Ami</b> Check Pulse VF/VT?	VF/VT? <b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse													
<b>Shock</b> CPR 2-min <b>Epi</b> Check Pulse	VF/VT? <b>Shock</b> CPR 2-min <b>Nothing</b> Check Pulse													
<b>H's</b>	Hypovolemia Hypoxia Hydrogen ion (acidosis) Hypo/Hyperkalemia Hypothermia													
<b>T's</b>	Tension Pneumothorax Tamponade (cardiac) Toxins Thrombosis (pulmonary) Thrombosis (coronary)													
Rate 100 - 120 Depth 2" Recoil No excessive ventilation Switch every 2 min MAP = 60														